




# VISIO-TEK™

Wireless digital microscope

THE PORTABLE  
DIGITAL MICROSCOPE  
WITH INTEGRATED WI-FI

## **RCKX:** Adjustable Vertical Rack for Microscope

The RCKX is a rigid vertical desktop rack that is designed to be used with all Visio-tek digital microscopes. It provides about 20cm of vertical working distance from the object.



Dimensions: (W)12.7cm x (L)20.7cm x (H)30.5cm

## **APPLICATIONS**

- *Industrial and Electronic Inspection/QC*
  - *Education*
  - *Research and Scientific Analysis*
  - *Forensics*
  - *Print inspection*
  - *Textile inspection*
  - *Trichology*
  - *Entomology*
  - *Botany*
  - *Mineralogy*
  - *Collectibles*
  - *Jewelry and Watch Repair*
  - *Reading assistance*
- And more ...*



[www.visio-tek.net](http://www.visio-tek.net) For more information: [info@visio-tek.net](mailto:info@visio-tek.net)

## TECHNICAL DATA SHEET