



# VISIO-TEK™

Wireless digital microscope

THE PORTABLE  
DIGITAL MICROSCOPE  
WITH INTEGRATED WI-FI

## **RCKBL:** BackLight Rack for Microscope

RCKBL is suitable for the complete line of Visio-tek digital microscopes.



- Up and down knobs are used to adjust the distance between microscope and object.
- With the removable X&Y crossed, you can slightly adjust the position of the object by rotating the knobs.
- Through the brightness control knob, you can choose to turn on or turn off the Led backlight source and adjust the brightness to observe the object much clear.

## **APPLICATIONS**

- *Industrial and Electronic Inspection/QC*
  - *Education*
  - *Research and Scientific Analysis*
  - *Forensics*
  - *Print inspection*
  - *Textile inspection*
  - *Trichology*
  - *Entomology*
  - *Botany*
  - *Mineralogy*
  - *Collectibles*
  - *Jewelry and Watch Repair*
  - *Reading assistance*
- And more ...*

