



# VISIO-TEK™

Wireless digital microscope

THE PORTABLE  
DIGITAL MICROSCOPE  
WITH INTEGRATED WI-FI

## VISIO-SCOPE5: Kit of VISIO-EYE Wi-Fi camera, VISIO-L04 Lens & Adapters and 3 Spacer Rings

VISIO-SCOPE5 is a complete solution with a digital C-Mount Wi-Fi camera (VISIO-EYE) that allows you to wirelessly digitalize any traditional optical device in seconds.

Visio-scope5 is compatible with all optical microscopes. Digitalize images and take measurements and photos from your standard optical device using your tablet, smartphone or PC with dedicated Apps and Software. The images can be displayed on different devices simultaneously.



## FEATURES

- Photo & video Wi-Fi mode resolution: 1.3MPx -1280x1024, 640x480, 320x240
- Photo & video USB mode resolution: 2592x1944 px, 2048x1536 px, 1600x1200 px, 1280x960 px, 1280x1024 px, 640x480 px, 320x240 px.
- Frame rate: up to 30fps
- IEEE 802.11 b/g/n standards compliant
- Wireless transmission distance: up to 5m
- Frequency range: 2.400GHz ~ 2.4835GHz
- Li-ion battery: continuous working time - appx. 2 hours  
full charging time - appx. 2 hours
- Camera sensor: CMOS 1/4"
- Lens: high definition microscopy lens
- Dimension: 36mm (∅) x 75mm (length)
- PC interface: USB 2.0
- Mini USB cable: 1.2m
- Power consumption: max 2.5W

## APPLICATIONS

- Industrial and Electronic Inspection/QC
  - Education
  - Research and Scientific Analysis
  - Forensics
  - Print inspection
  - Textile inspection
  - Trichology
  - Entomology
  - Botany
  - Mineralogy
  - Collectibles
  - Jewelry and Watch Repair
  - Reading assistance
- And more ...



www.visio-tek.net For more information: info@visio-tek.net

## TECHNICAL DATA SHEET