



VISIO-TEK™

Wireless digital microscope

THE PORTABLE
DIGITAL MICROSCOPE
WITH INTEGRATED WI-FI

VISIO-EYE5: Digital Wi-Fi camera (C-mount) for Optical Instruments & Devices + Spacer Ring

VISIO-EYE5 can be mounted without additional Lens or Adapters on any trinocular optical microscope with separate camera C-Mount port.



FEATURES

- Photo & video Wi-Fi mode resolution: 1.3MPx -1280x1024, 640x480, 320x240
- Photo & video USB mode resolution: 2592x1944 px, 2048x1536 px, 1600x1200 px, 1280x960 px, 1280x1024 px, 640x480 px, 320x240 px.
- Frame rate: up to 30fps
- IEEE 802.11 b/g/n standards compliant
- Wireless transmission distance: up to 5m
- Frequency range: 2.400GHz ~ 2.4835GHz
- Li-ion battery: continuous working time - appx. 2 hours
full charging time - appx. 2 hours
- Camera sensor: CMOS 1/4"
- Lens: high definition microscopy lens
- Dimension: 36mm (Ø) x 75mm (length)
- PC interface: USB 2.0
- Mini USB cable: 1.2m
- Power consumption: max 2.5W

APPLICATIONS

- Industrial and Electronic Inspection/QC
 - Education
 - Research and Scientific Analysis
 - Forensics
 - Print inspection
 - Textile inspection
 - Trichology
 - Entomology
 - Botany
 - Mineralogy
 - Collectibles
 - Jewelry and Watch Repair
 - Reading assistance
- And more ...



www.visio-tek.net For more information: info@visio-tek.net

TECHNICAL DATA SHEET