




VISIO-TEK™

Wireless digital microscope


THE PORTABLE
DIGITAL MICROSCOPE
WITH INTEGRATED WI-FI

VISIO-MOT: Digital Wi-Fi Otoscope, 1.3Mpx, 20x, 4 white LEDs

VISIO-MOT is a digital Otoscope that is light weight and offers portable flexibility in the examination area. Use with Visio-tek disposable speculum.



Ear Exam
Mild erythema



Ear Exam
Otitis

(FDA LISTED)

FEATURES

- Photo & video Wi-Fi mode resolution: 1.3MPx -1280x1024, 640x480, 320x240
- Magnification: 20x
- Object Distance: 1~160mm
- Frame rate: up to 30fps
- IEEE 802.11 b/g/n standards compliant
- Wireless transmission distance: up to 5m
- Frequency range: 2.400GHz ~ 2.4835GHz
- Li-ion battery: continuous working time - appx. 2 hours
full charging time - appx. 2 hours
- Camera sensor: CMOS 1/4"
- Focus range: 8mm ~ 60mm
- Lens: high definition microscopy lens
- Dimension: 36mm (Ø) x 142mm (length)
- PC interface: USB 2.0
- Mini USB cable: 1.2m
- Power consumption: max 2.5W
- Light source: 4 adjustable white LEDs

APPLICATIONS

- ENT
- Emergency Room
- Research & Scientific Analysis
- Telemedicine
- Telehealth
- Education



www.visio-tek.net For more information: info@visio-tek.net

TECHNICAL DATA SHEET